

Cleanroom Robotic Logistics


JÖRN SCHMID

GOLDFUSS ENGINEERING GMBH

RAYA 2025 Finalist Event



The Challenge for Cleanroom Robotic Logistics



- Cleanroom areas grade D to A
- **▶** Product logistics




Ready-To-Use Tubs



Connecting **Pharmaceutical** Knowledge

The Solution for Cleanroom Robotic Logistics



- Infeed airlock: Manual removal of outer packaging; RFID labelling of tubs (Class C)
- 2 Robot 1: Transfers tubs to logistics section (Class C/robot)
- **3 Robot 1:** Stores tubs in 4-unit groups after ID check (Class C)
- 4 Robot 1: Sends tubs through UV-C tunnel (airlock between Classes C, B, A)
- **Solution** Specification (Specification) **Solution** Specification (Specification) **Solution** (Specific
- **6** Conveyor: Transfers to filling lines (Class A)
- **7 Transfer:** Filled containers from Class A to B to Class C or D (1st floor)
- **8** Robot 1 Feeds output airlock with finished products.



Connecting Pharmaceutical Knowledge ispe.org

Video



Available after ISPE event!





Connecting Pharmaceutical Knowledge ispe.org 4

Features of the Cleanroom Robotic Logistics

Double functionality of both robots



Robot 1 uses a four-position tool with disinfectant for automated rack cleaning, increasing efficiency and safety without ladders or telescopic tools.



Robot 2 manages cyclic placement, opening, and replacement of sedimentation plates for cleanroom grade A (RABS), with plate magazines prefilled during line setup.



Connecting Pharmaceutical Knowledge ispe.org

Features of the Cleanroom Robotic Logistics



1. Industry Attractiveness

Fully automated, lights-out logistics with integrated buffering, decontamination, and sterile transfer. Stäubli robots deliver seamless aseptic processing at 2 tubs per minute.



2. Economic Potential

Saves 3 FTEs per shift (€180k/year) and doubles throughput. Smaller footprint and modular reusability deliver ROI < 5 years.



3. Pharma Applicability

Supports Fill & Finish with RTU containers (syringes, vials, carpules). Ideal for CMOs and pharma with small to medium batches and frequent changeovers.



4. Flexibility

Handles all RTU tub formats with random-access buffering and multi-line integration. RFID and format-free grippers minimize cleanroom interventions.



5. Ease of Integration

Applicable to standardized tubs without customization (GMP qualification required). Six validated syringe variants ensure broad transferability.



Connecting Pharmaceutical Knowledge ispe.or

FINAL MESSAGE



- Protecting Europe. Creating the Future.
- Together Unique. Together Strong.
- Pioneering Pharma-Robotics Worldwide.
- Solve Skills Shortage with Lights-Out Factories.

This project is just an example – many more will follow.

Challenge the machine builders!

